



POINT- OF-CARE

osom[®]
hCG Family Rapid Tests

OSOM[®] hCG Tests

FOR THE QUALITATIVE DETECTION OF HUMAN CHORIONIC GONADOTROPIN (hCG)

Features and Benefits

- **Accurate** – No biotin interference
- **Low detection limits** – 10 mIU/mL (serum) and 20 mIU/mL (urine)*
- **Fast** – 3 minutes* (urine) and 5 minutes (serum)
- **Easy** – Only one step. CLIA waived for Urine
- **Convenient** – Transfer pipette included with each test
- **Exceptional Customer Service** – Online training modules and technical support staffed by certified medical technologists/professionals

Performance Characteristics

| | Sample Type | Limit of Detection | Time to Result |
|--------------------------------|-------------|--------------------|----------------|
| hCG Urine dipstick (cat # 101) | Urine** | 25 mIU/mL | 5 minutes |
| hCG Urine Card (cat # 102) | Urine** | 20 mIU/mL | 3 minutes |
| hCG Combo Card (cat # 124) | Serum | 10 mIU/mL | 5 minutes |
| | Urine** | 20 mIU/mL | 3 minutes |

*cat # 102 and cat #124

**CLIA Waived for urine



OSOM® hCG Tests

SIMPLE TO USE / RESULTS IN MINUTES

Ordering Information

| | Configuration | Catalog Number | Storage | CPT Code |
|--------------------|----------------------|----------------|--------------|--------------------------------|
| hCG URINE DIPSTICK | 50 tests | 101 | Room Temp | Urine: 81025QW Serum: 84703 |
| hCG URINE CARD | 25 tests | 102 | Room Temp | |
| hCG COMBO CARD | 25 tests | 124 | Room Temp | |
| hCG URINE CONTROL | 2 X 10mL (Pos & Neg) | 134 | Refrigerated | |
| hCG SERUM CONTROL | 2 X 5mL (Pos & Neg) | 138 | Refrigerated | |

To learn more about all of our OSOM Tests scan the code or go to osomtests.com



THE AMERICAS

SEKISUI Diagnostics, LLC
One Wall Street
Burlington, MA 01803
Phone: 800 332 1042
Fax: 800 762 6311

Email: questions@sekisui-dx.com



SEKISUI
DIAGNOSTICS

Because every result matters™

sekisuidiagnostics.com