

# CURIE

## Versatility that Delivers



### Efficient Curing Space

- Perfect footprint for dental offices and labs.



### Future-proof

- Easy resin library update mechanism with SD card reader.



### Validated Workflows

- Over 150 preloaded resin profiles available at the touch of a button.



### Premium Performance

- 360-degree light distribution covering 365, 385 and 405 nm achieving superior curing results.



## CURIE Technical Specifications



### Curing Time Cycles

1 sec. to 20 min.



### Dimensions

7.9 (W) x 8.7 (D) x 6.3 (H) in



### Weight

6.6 lbs



### Chamber Dimensions

5.5 (W) x 5.5 (D) x 3.1 (H) in



### UV Wavelength

365, 385, 405 nm



### Technology

UV LED Array



### Power Consumption



75 W



### Power Input

100-240VAC, 50-60Hz

Item Number  
**CURIE - KT200200101**



## PACKAGE CONTENTS

- 1 x Glass Panel Box
- 1 x Power Cable
- 1 x Ethernet Cable
- 1 x Wi-Fi Dongle



Scan the QR Code to  
Schedule a Virtual Demo Today!