



SOL

Compatible resins by brand

List of FDA Approved Resin

June 2024

Available Build Platforms



Large Build Platform **L**

Printing Area: 128mm x 80mm
5.1 x 3.1 inches



Medium Build Platform **M**

70mm x 70mm
2.7 x 2.7 inches



Small Build Platform **S**

40mm x 20mm
1.6 x 0.8 inches

CURO

Resin Type:

Available Build
Platforms:



Cast

CURO Cast
CURO Pro Cast

L M S

Gingiva

CURO Gingiva

L M S

Guide

CURO Guide

L M

IBT

CURO I-B-T

L M

Temp

CURO Crown A1
CURO Crown A2
CURO Crown A3.5
CURO Crown B1
CURO Crown B3

L M S

Model

CURO Model
CURO Element Model Beige
CURO Element Model White
CURO Ortho Model

L M

Splint

CURO Splint
CURO Pro Splint

L M

Tray

CURO Tray

L M



Resin Type:

Available Build Platform:

Model

Argen Model Beige
Argen Model Light Grey



Tray

Argen Tray



Resin Type:

Available Build Platform:


Guide

Argen Guide



Splint

Argen Flex





Resin Type:

Available Build Platform:

Model

Model HP

L

Resin Type:

Available Build Platform:

Perm

VarseoSmile Crown plus

VarseoSmile TriniQ

L M S

Temp

VarseoSmile Temp

VarseoSmile Teeth

L M S

L M



Available Build Platform:

Cast

Dental Cast



Guide

Surgical Guide



Model

Model Standard
Model Basic



Temp

C&B



Available Build Platform:

Cast

3Delta Cast P
3Delta Cast M




Gingiva

3Delta Gingiva Mask



Model

3Delta Model 320
3Delta Model Grey



Resin Type:

Available Build Platform:

Resin Type:

Available Build Platform:

Denture

Dentca Denture Base II

L M

Temp

- Dentca Crown and Bridge A1
- Dentca Crown and Bridge A2
- Dentca Crown and Bridge A3
- Dentca Crown and Bridge A3.5
- Dentca Crown and Bridge B1
- Dentca Crown and Bridge B2

L M S


Try in

Dentca Try-In

L M

Dentca Denture Teeth

L M



**Available Build
Platform:**

Resin Type:

Cast

Optiprint Cast
Optiprint Zero
Optiprint Cast NF

L M S

Gingiva

Optiprint Gingiva

L M S

Guide

Optiprint Guide

L M



IBT

Optiprint I-B-T

L M

Model

Optiprint Model Golden Brown
Optiprint Model Grey
Optiprint Model Ivory
Optiprint OrthoModel

L M

Perm

Optiprint Lumina A1
Optiprint Lumina A2
Optiprint Lumina A3.5
Optiprint Lumina B3

L M S

Splint

Optiprint Clara
Optiprint Splint

L M

Temp

Optiprint Temp A2
Optiprint Temp A3
Optiprint Temp B2

L M S

Tray

Optiprint Tray

L M

Other

Optiprint Match
Optiprint Mock Up

L



Resin Type:

Model

- Luxaprint Model BGE
- Luxaprint Model GRY
- Luxaprint Model Transparent

Tray

- Luxaprint Tray

Gingiva

- Luxaprint Gingiva

Available Build Platform:



Resin Type:


Model

- FotoDent model2 beige
- FotoDent model2 white

Guide

- FotoDent guide

Available Build Platform:



Resin Type:

Available Build Platform:

Resin Type:

Available Build Platform:

Model



- KeyModel Ultra Ivory
- KeyModel Ultra Light Grey
- Keystone KeyModel
- Keystone KeyOrthoModel

Gingiva



KeyMask

Guide



KeyGuide

IBT



KeyOrtho IBT

Splint



- KeySplint Soft
- KeySplint Hard
- KeySplint Soft Clear
- KeySplint Hard Clear

Tray



KeyTray

Try in



Keystone Try-In





Resin Type:

Available Build Platform:

Available Build Platform:

Available Build Platform:

Cast



Mack4D Cast

Mack4D Tiko G

Denture



Mack4D Denture Dark Pink

Mack4D Denture Light Pink

Guide



Mack4D Guide

Model



Mack4D Model Golden Brown

Gingiva



Mack4D Gingiva

IBT



Mack4D I-B-T

Hearing



Hearing Cast

Mack 4D Audio-mold

Other



Mack4D Match

Mack4D Mock Up

Splint



Mack4D Splint

Temp



Mack4D Temp A1

Mack4D Temp A2

Mack4D Temp A3

Tray



Mack4D Tray





Available Build Platform:

Resin Type:

Cast

NextDent Cast



Denture

NextDent Denture 3D+



Guide

NextDent SG



Gingiva

NextDent Gingiva Mask



IBT

NextDent Ortho IBT



Available Build Platform:

Model

NextDent Model 2.0



Splint

NextDent Ortho Clear

NextDent Ortho Flex

NextDent Ortho Rigid



Temp

NextDent CMFH



Tray

NextDent Tray



Try In

NextDent Try In





Resin Type:

| Denture

Rodin Denture

| Guide

Rodin Surgical Guide

| IBT

Rodin Ortho IBT

| Model

Rodin Model

Available Build Platform:



Rodin Titan



Resin Type:

| Perm

Rodin Sculpture
Rodin Sculpture 2.0
Titan

| Splint

Rodin Splint

| Tray

Rodin Bite Tray

| Try In

Rodin Try-In

Available Build Platform:





Resin Type:

■ Perm

Saremco print CROWNTEC

Available Build
Platform:

L M S



Resin Type:**Available Build
Platform:****| Cast**

P pro Cast P
P pro Cast M

**| Gingiva**

P pro Gingiva Mask

**| Model**

P pro Master Model Dark Beige
P pro Master Model Grey





Resin Type:

Available Build Platform:

Resin Type:

Available Build Platform:

Cast

VOCO V-Print Cast



Denture

VOCO V-Print Dentbase



Guide

VOCO V-Print SG



Splint

VOCO V-Print Splint comfort
VOCO V-Print Splint



Model

VOCO V-Print Model Beige
VOCO V-Print Model Fast
VOCO V-Print Model 2.0



Temp

VOCO V-Print c&b temp



Tray

VOCO V-Print Tray



Try in

VOCO V-Print Try-In





Resin Type:

Cast

VeriCast

Gingiva

VeriGUM

Guide

VeriGuide

Tray

VeriModel Tray

Available Build Platform:



Resin Type:

Model

VeriModel Golden Brown

VeriModel Gray


VeriModel Ivory

VeriModel White

Try In

VeriDent Try-In

Available Build Platform:





SOL

Compatible resins by brand

List of Health Canada Approved Resin

June 2024

Available Build Platforms



Large Build Platform **L**

Printing Area: 128mm x 80mm
5.1 x 3.1 inches



Medium Build Platform **M**

70mm x 70mm
2.7 x 2.7 inches



Small Build Platform **S**

40mm x 20mm
1.6 x 0.8 inches



Resin Type:

Available Build
Platforms:

Model

CURO Model


CURO Element Model Beige

CURO Element Model White



Gingiva

CURO Gingiva





Resin Type:

Available Build Platform:

Model

Argen Model Beige
Argen Model Light Grey



Resin Type:

Available Build Platform:

Guide

Argen Guide



Tray

Argen Tray



Splint

Argen Flex





Resin Type:

**Available Build
Platform:**

Temp

VarseoSmile Temp

VarseoSmile Teeth

L M S



**Available Build
Platform:**

Cast

3Delta Cast P
3Delta Cast M




Gingiva

3Delta Gingiva Mask



Model

3Delta Model 320
3Delta Model Grey



Resin Type:

Available Build Platform:

Resin Type:

Available Build Platform:

Denture

Dentca Denture Base II

L M

Temp

- Dentca Crown and Bridge A1
- Dentca Crown and Bridge A2
- Dentca Crown and Bridge A3
- Dentca Crown and Bridge A3.5
- Dentca Crown and Bridge B1
- Dentca Crown and Bridge B2

L M S


Try in

Dentca Try-In

L M

Dentca Denture Teeth

L M



Resin Type:

Model

KeyModel Ultra Ivory
KeyModel Ultra Light Grey
Keystone KeyModel
Keystone KeyOrthoModel

Splint

KeySplint Soft
KeySplint Hard
KeySplint Soft Clear
KeySplint Hard Clear

Available Build Platform:



Resin Type:

Gingiva

KeyMask

Guide

KeyGuide

IBT

KeyOrtho IBT


Tray

KeyTray

Try in

Keystone Try-In

Available Build Platform:





Available Build Platform:

Resin Type:

Cast

NextDent Cast



Denture

NextDent Denture 3D+



Guide

NextDent SG



Gingiva

NextDent Gingiva Mask



IBT

NextDent Ortho IBT



Available Build Platform:

Model

NextDent Model 2.0



Splint

NextDent Ortho Clear

NextDent Ortho Flex

NextDent Ortho Rigid



Temp

NextDent CMFH



Tray

NextDent Tray



Try In

NextDent Try In



Resin Type:

| Denture

Rodin Denture

| Guide

Rodin Surgical Guide

| IBT

Rodin Ortho IBT

| Model

Rodin Model

Available Build Platform:



Resin Type:

| Splint

Rodin Splint

| Tray

Rodin Bite Tray

| Try In

Rodin Try-In

Available Build Platform:



Resin Type:**Available Build
Platform:****| Cast**

P pro Cast P
P pro Cast M

**| Gingiva**

P pro Gingiva Mask

**| Model**

P pro Master Model Dark Beige
P pro Master Model Grey





Resin Type:

Available Build Platform:

Resin Type:

Available Build Platform:

Cast

VOCO V-Print Cast



Denture

VOCO V-Print Dentbase



Guide

VOCO V-Print SG



Splint

VOCO V-Print Splint comfort
VOCO V-Print Splint



Model

VOCO V-Print Model Beige
VOCO V-Print Model Fast
VOCO V-Print Model 2.0



Temp

VOCO V-Print c&b temp



Tray

VOCO V-Print Tray



Try in

VOCO V-Print Try-In

