


lightforce®  
therapy lasers



## *i*nfluence® Technology

*i*nfluence treatments. *i*nfluence patients. *i*nfluence outcomes.

**/// EXPi & /// FXi**

## CUSTOM



**influence** technology is patient specific.

- Condition
- Body Type
- Skin Type
- Acuity

SELECT PATIENT CHARACTERISTICS

Height	Body Type	Color
Short <input checked="" type="radio"/>	Small <input type="radio"/>	Light (I, II) <input type="radio"/>
Medium <input type="radio"/>	Medium <input checked="" type="radio"/>	Medium (III, IV) <input checked="" type="radio"/>
Tall <input type="radio"/>	Large <input type="radio"/>	Dark (V, VI) <input type="radio"/>

Fitzpatrick Scale

< Previous Next >

## What is **influence** technology?

The newest innovation from the leaders in laser therapy. Harnessing software and hardware advancements, LightForce® brings this "next level" technology to new heights and delivers outcomes like never before. Dose and delivery are prime factors and can vary based on patient specific characteristics. The NEW **i**-series provides utility where and when you need it.

# Get in contact and

## CONTACT



**influence** technology is different.

- Patented Delivery
- Reach Target Tissue
- Unmatched Versatility
- Manual Application

## Versatile delivery

**influence** outcomes with the ability to deliver treatments in both contact and non-contact delivery methods.

**EMPOWER**  
Delivery System™





CONSISTENT



*influence* technology is repeatable.

- Save Settings
- Recall Favorites
- Reliable Results
- Easy to Use

## *influence*® Technology Deep Tissue Laser Therapy™

*influence* treatments. *influence* patients. *influence* outcomes.

and achieve results

### ▶ Benefits of contact application



**Compression** - Reach targeted tissue by removing superficial absorbers



**Collimation** - Tighter beam minimizes energy loss



**Reflection** - Significantly reduced in contact mode



**Soft Tissue** - Manually work tissue while delivering energy



**Refractive Index** - Minimize light loss due to similar index of skin and lens composition

The time with your patients is valuable

*influence* technology will maximize your treatments and provide the consistent results you expect with your patients. Treat with confidence knowing you are making the most out of your treatment time.



5-7 min

## The LightForce<sup>®</sup> *i* Series

### *EXPi*

25 Watts



- Fast and efficient treatments
- High power industry standard

**Laser Type:** Class IV, Solid State

**Wavelengths:** 980/810 nm

**Operating Modes:** CW or Pulsed

**Dimensions:** 15.5"x9.75"x10.5"

**Weight:** 17 lbs

**Warranty:** 24 months

### *FXi*

15 Watts



- Functional mobility, and battery operation
- Price performance leader

**Laser Type:** Class IV, Solid State

**Wavelengths:** 980/810 nm

**Operating Modes:** CW or Pulsed

**Dimensions:** 15"x11"x 8"

**Weight:** 7 lbs

**Warranty:** 24 months

### EDS

Empower  
Delivery  
System



- 5 application-specific treatment attachments including the patented massage ball applicator
- Premium, double sheathed, rubber coated fiber optic cable
- Integrated finger switch
- 24 month warranty

Contact LightForce<sup>®</sup> Therapy Lasers

@ LightForceLasers.com

☎ 877-627-3858

✉ 250 Corporate Blvd, Suite B, Newark, DE 19702

**lightforce**<sup>®</sup>  
therapy lasers