



Upgrade Your Point-of-Care Hemoglobin Testing with the HemoCue® Hb 801 System

	HemoCue® Hb 201+ System	HemoCue® Hb 801 System
Speed	Results can take up to 60 seconds	Results in less than 1 second
Microcuvettes	Reagent microcuvettes	Reagent-free microcuvettes for improved open vial stability
Microcuvette shelf life	18 months	30 months
Microcuvette stability	3 months for open vials	30 months for open vials *From date of manufacturing
Display	Small display, error codes for maintenance	Large display with intuitive on-screen symbols
Connectivity	None	Bluetooth and middleware option *Not approved for hospital use




CLIA WAIVED



CLIA WAIVED


Accuracy Always Matters

Hemoglobin testing is often done as a part of a routine doctor's visit or when a patient presents symptoms of anemia.¹ With a small drop of blood (~10 µL), the HemoCue® Hb 801 System determines an individual's hemoglobin value in less than one second. Calibrated against the International Council for Standardization in Hematology (ICSH) reference method, the HemoCue® Hb 801 System provides lab-quality results in CLIA-waived settings.



Each point on the diagram represents an actual hemoglobin test result. The closer to the line, the more closely related the result is to the "gold standard" reference method.²

HemoCue OnCue™ Support



One-Touch Resolution

Get answers from our live customer service and technical support teams located at HemoCue® America in California.

In-Field or Virtual Service

Get extensive support from your HemoCue® representative virtually or right in your office.

Online Education

Easily stay up to date with comprehensive eLearning modules, competency reviews, and more.

Committed to Accuracy Since Day One

Forty years ago, we set out to do the impossible — bring lab-quality hemoglobin testing to the point of care. Today, **more healthcare professionals still choose HemoCue** systems for point-of-care hemoglobin testing than systems from any other brand.³

Product Order Number

121920

HemoCue® Hb 801 Analyzer

111902

HemoCue® Hb 801 Microcuvettes

Make the HemoCue® Hb 801 System part of your workflow.

Contact your HemoCue or preferred distribution representative.

Learn more at hemocue.us



¹ <https://www.mayoclinic.org/tests-procedures/hemoglobin-test/about/pac-20385075>

² HemoCue® Hb 801 System V&V Graphs, Figure 2; HCAB:000020941[2]

³ Sanchis-Gomar et al. Hemoglobin Point-of-Care Testing: The HemoCue System