



COPY DENTURES

Copy dentures can be fabricated for an upper or lower denture. The main requirement is that the scans be a full 360 degrees of coverage and no tears or holes. The appliances are ideal for those patients who insist on a new denture just like their old denture. There are limitations:

- We cannot affect or change the VDO
- We are limited to giving the patient exactly what they currently have with new acrylic and teeth
- The anatomy of the teeth remain the same
- In most cases we may thicken the area around the necks of the teeth in order to create tooth pockets for assembly

Upper Sample Scan



Lower Sample Scan

